

Assistance with settings for the Quantisal Collector Training:

If the collector training window pops up but there are no navigation buttons (ie next/back) at the bottom of the screen:

- 1) Make sure the Zoom Setting is at 100% (or less): At the Internet Explorer menu bar go to
 - View
 - Zoom
 - Select 100%

If this still doesn't work adjust the resolution on your computer:

- Right click anywhere on her desktop and select 'Properties'. (This is usually the last option in the menu that pops up.)
- Click the 'Settings' tab.
- Drag the slider that is on the lower left part of the window to the desired resolution. In this case, 800x600.